Instructions for use of the procedure pack

Description

A procedure pack designed for isolation and separation of platelet-rich plasma (PRP) from the patient's blood, and then its simultaneous intra-articular injection with sodium hyaluronate.

Monoshot HA with pH 6.8-7.4 contains high molecular weight sodium hyaluronate produced by bacterial fermentation, non-animal origin. It is sterilized and non-pyrogenic. Monoshot HA is intended for intra-articular injection to increase the viscosity of the synovial fluid.

Monoshot PRP medical device is used for the isolation and separation of platelet-rich plasma (PRP) from patient's blood. Isolated autologous PRP is intended for local injections to improve the regenerative properties of cells and tissues by providing a large number of growth factors to specific areas. Simultaneous injection of PRP and sodium hyaluronate supports regenerative processes, and their synergistic action may result in reduction of the inflammatory response, thus leading to pain reduction and improvement of functionality of the affected joint. Use of the set designed in this way reduces the duration of the operation procedure and leads to lower traumatization of the patient.

Procedure pack consists of the following medical devices:

- ① Monoshot PRP, *Biovico* (separation device)
- ② Venofix® A, 21G, *Nipro Corporation* (butterfly needle)
- ③ Microlance™ 3, 21G x 1½" – 2 pcs., *BD* (needle)
- ④ dicoNEX Single use syringe, 3-part 20 ml, *Zarys*
- ⑤ dicoNEX Single use syringe, 3-part 10 ml, *Zarys*
- ⑥ Microlance™ 3, 19G x 1½", *BD* (needle)
- ⑦ Monoshot HA, *Biovico* (3.8 ml of 2.5% sodium hyaluronate in pre-filled syringe)

Indications

Pain and decreased articular mobility associated with traumatic or degenerative changes including osteoarthritis, where there is a need to improve the regenerative properties at the specific areas. Due to the administration of high-concentrated sodium hyaluronate in combination with PRP in total volume of 4.8 ml, the procedure pack is dedicated mostly for knee joints.

Contraindications

The procedure pack should not be administered to patients:

- with known hypersensitivity to any of the components of Monoshot HA;
- with septic arthritis;
- with dermatosis at the injection site;
- who are using anticoagulants.

Children, pregnant or nursing women should not be treated as there are no clinical data available on therapy in those groups of patients.

Precautions

To be used exclusively by qualified, trained personnel. Before using the procedure pack, please read the instructions for use. The procedure pack consists of disposable medical devices. Do not reuse. Do not use any of medical devices in case of damage of its direct package.

Medical devices are sterile only during their shelf life indicated on the package.

Store the procedure pack at temperature between +5°C and +25°C. Protect from direct sunlight, high temperature and high humidity. Keep out of the reach of children.

Wear personal protective equipment when using the procedure pack. Ensure that no one is exposed to direct contact with collected blood. Place all used pack parts in a medical waste container and utilize them.

Do not concomitantly use disinfectants containing quaternary ammonium salts for skin preparation as hyaluronan can precipitate in their presence.

Local symptoms such as pain, redness and swelling may occur at the injection site. These effects can be reduced by cooling the site with cooling bags for 5 to 10 minutes after the injection. Simultaneous oral intake of analgesics and anti-inflammatory drugs (NSAIDs) may be beneficial for pain relief.

NOTE: Use only freshly drawn blood for plasma isolation. PRP must be administered immediately after isolation.

Preparation before use

1. Remove Monoshot PRP ① from the package, position it vertically and mark with the patient's code or name to avoid confusions.
2. Using the needle ③, withdraw 1.5 ml of anticoagulant into 20 ml syringe ④.
3. Using butterfly needle ②, collect the patient's blood into syringe containing anticoagulant. Fill the syringe up to 15 ml.
4. Mix gently by inversion.
5. Mark syringe containing collected blood with the patient's code or name.

NOTE: To avoid contamination of the isolated PRP, perform the entire isolation procedure with caution to sterility.

Instructions for use

1. Check if the bottom screw cap of Monoshot PRP ① is tightened correctly. Transfer carefully collected blood with anticoagulant using the needle ③ into Monoshot PRP ① placed vertically on a flat surface, pushing the needle in the center of the top cap.
2. Open the centrifuge cover, place Monoshot PRP ① and properly prepared counterweights inside.
3. Centrifuge for **5 minutes at 1700 x g** (RCF).
4. Open the cover and carefully take out Monoshot PRP ① with centrifuged plasma. Hold Monoshot PRP ① upright, do not shake it in order to avoid mixing of separated fractions.
5. Replace the bottom screw cap of Monoshot PRP ① with controller's screw cap and then fasten the controller's bolt to Monoshot PRP ①.
6. Turn the bottom part of controller's bolt counterclockwise, until plasma reaches the last ring in the top part of Monoshot PRP ①.
7. Open the top cap of Monoshot PRP ①.

NOTE: Each time before connecting the syringe to Monoshot PRP ①, check if syringe plunger is moving properly.

8. Attach 10 ml syringe ⑤ to a place of top cap to collect platelet-poor plasma (PPP). Rotate a controller's bolt until red blood cells (RBC) fraction reaches the line marked as '1' on Monoshot PRP ①.
9. Take the pre-filled syringe Monoshot HA ⑦ out of the blister and remove Luer-lock adapter. Then, to withdraw platelet-rich plasma (PRP), connect the pre-filled syringe Monoshot HA ⑦ to the top cap and carefully turn the controller's bolt until the RBC fraction reaches the last ring in the top part of Monoshot PRP ①.
10. Disconnect the pre-filled syringe Monoshot HA ⑦ supplemented with the withdrawn PRP.
11. Mark the pre-filled syringe Monoshot HA ⑦ supplemented with PRP with the patient's code or name in order to avoid mistakes.
12. Disinfect the skin at the planned injection site. Screw a sterile needle ⑥ onto the pre-filled syringe ⑦. Remove any air bubbles before injection.

The therapy consists of one intra-articular injection. Effects of the therapy last for at least six months. The therapy can be repeated as needed.

Properties and mode of action

Sodium hyaluronate is responsible for the viscosity and elasticity of the synovial fluid. Because of its lubricating and shock-absorbing properties, it enables painless, physiological joint movement. Moreover, it also participates in the nutrition of the cartilage. Supplementation of synovial fluid by intra-articular injection of sodium hyaluronate, improves the viscoelastic properties of the synovial fluid. After one injection, joint mobility is improved and pain associated with degenerative changes is relieved.

Autologous platelet-rich plasma is used to improve cell regeneration and accelerate the healing process by providing large number of growth factors to specific areas. PRP supports mechanical properties of sodium hyaluronate by providing regenerative properties and enriching its composition with platelets derived growth factors. Nevertheless, the combination of the properties of sodium hyaluronate with the properties of autologous platelet-rich plasma is possible only when hyaluronic acid has a sufficiently high concentration, viscoelastic properties and molecular weight, while PRP contains the maximum concentration of platelets in the possible lowest volume.

Date of issue: 03/2023

Rev. 01/2023

Organization assembling procedure pack:

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